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The purpose of this document is to provide an advanced verification to the ASL 5000 lung simulator. The advanced verification provides additional testing and verification to the standard calibration document. This document includes verification of the Standard ASL5000 3L, 5L, 6L, and Preemie ASL attachment.

Prüfer	
Person in charge	
Unterschrift Signature	
Name	
Datum / Date	

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Changed firmware for TAI testing purpose then changed back to 4.7.99 before shipping

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Pressure Sensor Verification – Temperature drift uncertainty

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This process will verify the voltage output of the pressure sensor such that there is a minimum drift uncertainty during operation. The test will compare the amount of change to the pressure sensor voltage output at the sensor as well as in the embedded firmware program.

Testing utilizes onboard software called ad_test.exe and compares the drift to a NIST traceable digital multimeter.

This process is to record the relevant hardware and software versions that the customer is using such that the verification testing will be based on the correct versions.

Drucksensorverifizierung – Temperaturdriftunsicherheit Pressure Sensor Verification – Temperature drift uncertainty (2.1.1)									
Time (min.)	0	15	30	45	60	Max- Min	Toleranz tolerance (v)	Ergebnis Result	
PuTTY Terminal (V)							0,002	<input type="checkbox"/> PASS	<input type="checkbox"/> FAIL
Multi- meter (V)							0,002	<input type="checkbox"/> PASS	<input type="checkbox"/> FAIL
Eingesetzte Prüfmittel Test equipment									
Gerät / Device						ID	Cal-next		

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This process will verify the pressure and volume deviations within a simulation at the end of the expiration or the idle time before the next breath. The breath rate is set to 8BPM (or 7.5 s per breath). This verification analyzes the idle time of the last 1.5 s of the expiratory phase.

Drucksensorverifizierung – Statistik der Drucknullkalibrierung									
Pressure Sensor Verification – Pressure zero calibration statistic (2.1.2)									
Breath	1	2	3	4	5	Max-Min	Toleranz tolerance	Ergebnis Result	
Airway Pressure Offset [cmH ₂ O]							0,1 cmH ₂ O	<input type="checkbox"/> PASS	<input type="checkbox"/> FAIL
Volume Offset [ml]							10 mL	<input type="checkbox"/> PASS	<input type="checkbox"/> FAIL

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Volume Verification with Calibrated Syringe

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This process will verify the volume accuracy with three (3) NIST traceable syringes.

Approximately 5 strokes will be made with each syringe based on the type of testing (3L, 5L, 6L, or Preemie). The 100 mL and 10 mL syringes utilize an additional ~18 mL of compressible gas as the pressure release adapter.

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This process will verify the calculated compliance for a range of compliance values. This covers the 6L, 5L, 3L, and Premie ASL options. Compliance values of 50, 20, 3, 1, and 0.5 mL/cmH₂O will be verified based on the Standard (6L, 5L, 3L) or Premie configurations.

The 100 mL and 10 mL syringes utilize an additional ~18 mL of compressible gas as the pressure release adapter

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Compliance verification

(2.1.4)

	Standard ASL Configuration			Preemie ASL Attachment Configuration (Option)	
Syringe Discharge [mL]	1000	100	10	10	10
Reported Pressure [cmH ₂ O]				N/A	N/A
Set C [mL/cmH ₂ O] (incl. compensations, accounts for Compressible gas volume in "deadspace")					
Parasity C [mL/cmH ₂ O] (from additional compressible volume in the "stroke")	1	0.1	0.01	0.01	0.01
Volume "lost" [mL] (due to parasitary compliance)(reported P*Parasitic C)				N/A	N/A
Adjusted Volume Discharge [mL] (Piston displacement)(Syringe Discharge - Volume "lost")				N/A	N/A
C Actual [mL/cmH ₂ O] (Adjusted Volume / Reported P)				N/A	N/A
Error [%] from Set C				N/A	N/A
Tolerance	5% \triangle 0,2 mL/cmH ₂ O	5% \triangle 0,2mL/cmH ₂ O	5% \triangle 0,2 mL/cmH ₂ O	N/A	N/A
Result	<input type="checkbox"/> PASS <input type="checkbox"/> FAIL	<input type="checkbox"/> PASS <input type="checkbox"/> FAIL	<input type="checkbox"/> PASS <input type="checkbox"/> FAIL	N/A	N/A

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This process will verify the response of the ASL5000 to a spontaneous muscle pressure. The goal is to assure the error between the calculated model and the measured piston movement is minimized.

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This process will verify the dynamic response of the ASL5000 to a flow-pump model with varying breath rates (5, 10, 12, 15, 20, 25).

Standard ASL Configuration (500mL)					
Standard ASL-Konfiguration (500mL) (2.1.6a)					
Breath Rate [BPM]	Gemessen measured	Abweichung deviation	Toleranz (%) tolerance	Ergebnis Result	
5			2% von 500 mL	<input type="checkbox"/> PASS	<input type="checkbox"/> FAIL
10			2% von 500 mL	<input type="checkbox"/> PASS	<input type="checkbox"/> FAIL
12			2% von 500 mL	<input type="checkbox"/> PASS	<input type="checkbox"/> FAIL
15			2% von 500 mL	<input type="checkbox"/> PASS	<input type="checkbox"/> FAIL
20			2% von 500 mL	<input type="checkbox"/> PASS	<input type="checkbox"/> FAIL
25			2% von 500 mL	<input type="checkbox"/> PASS	<input type="checkbox"/> FAIL

Premie ASL attachment Configuration (100mL)					
Premie ASL-Aufsatz Konfiguration (100mL) (2.1.6b)					
Breath Rate [BPM]	Gemessen measured	Abweichung deviation	Toleranz (%) tolerance	Ergebnis Result	
5			2% von 100 mL	<input type="checkbox"/> PASS	<input type="checkbox"/> FAIL
10			2% von 100 mL	<input type="checkbox"/> PASS	<input type="checkbox"/> FAIL
12			2% von 100 mL	<input type="checkbox"/> PASS	<input type="checkbox"/> FAIL
15			2% von 100 mL	<input type="checkbox"/> PASS	<input type="checkbox"/> FAIL
20			2% von 100 mL	<input type="checkbox"/> PASS	<input type="checkbox"/> FAIL
25			2% von 100 mL	<input type="checkbox"/> PASS	<input type="checkbox"/> FAIL

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This process verifies piston accuracy at high flow rates. These rates differ based on the ASL5000 configuration used in the testing (6L, 5L, 3L, and Premie).

Spitzenwert des Durchflusses Peak Flow Standard ASL Configuration (2.1.7a)					
	Gemessen measured	Abweichung deviation	Toleranz (%) tolerance	Ergebnis Result	
Model Peak Flow [L/min]			N/A	N/A	N/A
Piston Peak Flow [L/min]			N/A	N/A	N/A
Model Peak Volume [mL]			N/A	N/A	N/A
Piston Peak Volume [mL]			1 %	<input type="checkbox"/> PASS	<input type="checkbox"/> FAIL
Piston/Model Difference (Flow) [L/min]			1 %	<input type="checkbox"/> PASS	<input type="checkbox"/> FAIL

Spitzenwert des Durchflusses Peak Flow Premie ASL Attachment Configuration (2.1.7b)					
	Gemessen measured	Abweichung deviation	Toleranz (%) tolerance	Ergebnis Result	
Model Peak Flow [L/min]			N/A	N/A	N/A
Piston Peak Flow [L/min]			N/A	N/A	N/A
Model Peak Volume [mL]			N/A	N/A	N/A
Piston Peak Volume [mL]			1 %	<input type="checkbox"/> PASS	<input type="checkbox"/> FAIL
Piston/Model Difference (Flow) [L/min]			1 %	<input type="checkbox"/> PASS	<input type="checkbox"/> FAIL

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This process verifies the resistance settings on the ASL5000. The resistance verification is based on a step function input and utilization of the RC time constant, where RC is equal to 63.2% of the maximum tidal volume.

Resistance Verification Widerstandsprüfung (2.1.8)					
Set Resistance [cmH ₂ O/L/s]	5	20	50	200	500
Compliance [mL/cmH ₂ O]	100	20	10	2.5	1
V _t [mL] (max volume from spreadsheet)					
V _t @ RC _t (V _t * 0.632)					
RC [ms] (time at V _t @RC)					
Calculated R [cmH ₂ O/L/s] (RC / Compliance)					
Abweichung (%) deviation					
Toleranz (%) Tolerance based on ASTM F-1100	12%	12%	12%	12%	12%
Ergebnis Result	<input type="checkbox"/> PASS	<input type="checkbox"/> PASS	<input type="checkbox"/> PASS	<input type="checkbox"/> PASS	<input type="checkbox"/> PASS
	<input type="checkbox"/> FAIL	<input type="checkbox"/> FAIL	<input type="checkbox"/> FAIL	<input type="checkbox"/> FAIL	<input type="checkbox"/> FAIL